



ERYBANK[®]	ANTI-H LECTIN
Application	For H antigen determination and subgrouping
Titre	~ 1:16 with O group cells
Class of Agglutinin	Phytohaemagglutinin
Source	Ulex europaeus seed extract
Specificity	Negative reacting with 'O _h ' phenotype
Reactivity	Reactivity pattern with different red cells, O > A ₂ > A ₂ B > B > A ₁ > A ₁ B
Standardisation	Follows AABB and FDA guidelines
Storage / Stability	2-8°C, 24 months

Presentation	Pack	Cat.No.
ANTI - H Lectin reagent	5 ml	10210005
	10 ml	10210010

