



<b>ERYBANK®</b>	<b>BOVINE SERUM ALBUMIN</b>
<b>Application</b>	Bovine Serum Albumin for serological applications
<b>Total Protein</b>	~ 22 % (g/dl)
<b>Conductivity</b>	6 - 8 x 10 <sup>-4</sup> Ohms/cm at 25+/- 1°C
<b>Source</b>	Selected Raw Bovine Serum
<b>Purity</b>	Globulin free, single band on gel electrophoresis
<b>Standardisation</b>	Follows AABB and FDA guidelines
<b>Storage / Stability</b>	2-8°C, 36 months

<b>Presentation</b>	<b>Pack</b>	<b>Cat.No.</b>
<b>Bovine Serum Albumin reagent</b>	5 ml	10230005
	10 x 5 ml	10230105
	10 ml	10230010
	6 x 10 ml	10230610

