



**ERYCARD™ ABO** **Blood Grouping Card Test for ABO Forward Grouping with Autocontrol**

<b>Application</b>	<ul style="list-style-type: none"> <li>• Bed side ABO testing of Patient and Donor's blood unit before transfusion</li> <li>• Preconfirmation test, especially at nursing homes</li> <li>• For releasing compatible blood for emergency transfusions in             <ul style="list-style-type: none"> <li>• Defence Camps</li> <li>• Ambulance</li> <li>• Casualty services</li> </ul> </li> </ul>
--------------------	---

<b>Principle</b>	Self performing, based on the principle of lateral flow guided by capillary action
------------------	--

<b>Sensitivity &amp; Specificity</b>	100% correlation with classical tube technique
--------------------------------------	--

<b>Reagents</b>	<p>Reagents are pre-dried at the appropriate sample pad</p> <p>Anti-A (IgM) monoclonal antibodies in sample well A</p> <p>Anti-B (IgM) monoclonal antibodies in sample well B</p> <p>Autocontrol without any antibodies (Ctrl)</p> <p>Reagent Buffer (Proprietary Buffer)</p>
-----------------	---

<b>Standardisation</b>	Follows FDA guidelines
------------------------	------------------------

<b>Test duration</b>	Less than five minutes
----------------------	------------------------

<b>Storage / Stability</b>	2-30°C, 24 months
----------------------------	-------------------

Presentation	Pack	Cat.No.

<b>ERYCARD™ ABO</b> kit contains	24 Cards	102840024
Individually pouched device		
Reagent buffer		
Sample collection loops		

