



<b>ERYCLONE<sup>®</sup></b>	<b>ANTI-A</b>
<b>Application</b>	ABO grouping reagents for slide and tube tests
<b>Titre</b>	≥ 1: 256
<b>Class of Antibody</b>	IgM, Murine Monoclonal
<b>Source</b>	Murine hybridoma
<b>Specificity</b>	ANTI-A - 100% to A <sub>1</sub> , A <sub>2</sub> and A <sub>x</sub> antigens
<b>Standardisation</b>	Follows AABB and FDA guidelines
<b>Storage / Stability</b>	2-8°C, 24 months

Presentation	Cat.No.				
	ABO grouping reagents	Pack	5 ml	10x5 ml	10 ml
ANTI-A		101620005	101620105	101620010	101620610