



ERYCLONE[®]	ANTI-A,B
Application	ABO grouping reagents for slide and tube tests
Titre	≥ 1: 256
Class of Antibody	IgM, Murine Monoclonal
Source	Murine hybridoma
Specificity	ANTI-A,B - 100% to A and B antigens
Standardisation	Follows AABB and FDA guidelines
Storage / Stability	2-8°C, 24 months

Presentation	Cat. No.				
	ABO grouping reagents	Pack	5 ml	10x5 ml	10 ml
	ANTI-A,B	101640005	101640105	101640010	101640610