



	ANTI-C+D+E
Application	Slide & tube test reagent for Rh genotyping
Titre	<ul style="list-style-type: none"> • ~ 1:32 with r¹r (dCe/dce) phenotype • ~ 1:32 with R₁r (Dce/dce) phenotype • ~ 1:32 with r''r (dcE/dce) phenotype
Class of Antibody	IgM + IgG for D antigen, IgM for C,E antigen
Source	Human cell line
Specificity	<ul style="list-style-type: none"> • 100% to C, D, E antigen
Standardisation	Follows AABB and FDA guidelines
Storage / Stability	2-8°C, 24 months

Presentation	Pack	Cat.No.
Anti C+D+E	5 ml	10174005

