



|                             |  |
|-----------------------------|--|
| <b>ERYCLONE<sup>®</sup></b> | <b>ANTI-D</b>  |
| <b>Application</b>          | Rh <sub>0</sub> (D) typing reagent for slide and modified tube tests |
| <b>Titre</b>                | ≥ 1: 32  |
| <b>Class of Antibody</b>    | IgG, Monoclonal  |
| <b>Source</b>               | Human B cell line, EBV transformed                                   |
| <b>Specificity</b>          | 100% to Rh <sub>0</sub> (D) antigen                                  |
| <b>Standardisation</b>      | Follows AABB and FDA guidelines                                      |
| <b>Storage / Stability</b>  | 2-8°C, 24 months   |

| <b>Presentation</b>                     | <b>Pack</b> | <b>Cat.No.</b> |
|---|-------------|----------------|
| <b>Rh<sub>0</sub>(D) typing reagent</b> | 5 ml        | 10140005       |
|   | 10 x 5 ml   | 10140105       |
|   | 10 ml       | 10140010       |
|   | 6 x 10 ml   | 10140610       |

