



ERYCLONE[®]	ANTI-D
Application	Saline reacting Rh ₀ (D) typing reagent for slide and tube tests
Titre	≥ 1: 256
Class of Antibody	IgM, Monoclonal
Source	Human B cell line, EBV transformed
Specificity	100% to Rh ₀ (D) antigen
Standardisation	Follows AABB and FDA guidelines
Storage / Stability	2-8°C, 24 months

Presentation	Pack	Cat.No.
Rh₀(D) typing reagent	5 ml	101650005
	10 x 5 ml	101650105
	10 ml	101650010
	6 x 10 ml	101650610

