



ERYCLONE[®]	ANTI-E
Application	Slide & tube test Reagent for Rh genotyping
Titre	<ul style="list-style-type: none"> • ~ 1:32 with r'r (dCe/dce) phenotype • ~ 1:32 with rr (dce/dce) phenotype • ~ 1:32 with r''r (dcE/dce) phenotype
Class of Antibody	IgM, Monoclonals
Source	Human cell line
Specificity	<ul style="list-style-type: none"> • ANTI-E 100% to E antigen
Standardisation	Follows AABB and FDA guidelines
Storage / Stability	2-8°C, 24 months

Presentation		Pack	Cat.No.
Rh genotyping reagents	ANTI-E	2 ml	10172002
		5 ml	10172005

