



ERYCLONE[®]	ANTI HUMAN GLOBULIN
Application	Polyspecific reagent for Direct and Indirect AHG (Coomb's) Test
Titre	<ul style="list-style-type: none"> • ANTI-IgG ~ 1: 256 • ANTI-C3d ~ 1:16
Class of Antibody	<ul style="list-style-type: none"> • ANTI-IgG, Purified goat antibodies • ANTI-C3d, IgM murine monoclonal
Source	Goat serum / Murine hybridoma
Specificity	Polyspecific with activity against human IgG, and complement components C3b and C3d
Standardisation	Follows AABB and FDA guidelines
Storage / Stability	2-8°C, 24 months

Presentation	Pack	Cat.No.
Polyspecific AHG reagent (Green)	2 ml	10180002
	5 ml	10180005
	10 ml	10180010
	10 x 5 ml	10180105
	6 x 10 ml	10180610

