



<b>ERYCLONE<sup>®</sup></b>	<b>ANTI-c</b>
<b>Application</b>	Slide & tube test Reagent for Rh genotyping
<b>Titre</b>	<ul style="list-style-type: none"> <li>• ~ 1:32 with r<sup>1</sup>r (dCe/dce) phenotype</li> <li>• ~ 1:32 with rr (dce/dce) phenotype</li> <li>• ~ 1:32 with r<sup>1</sup>r (dcE/dce) phenotype</li> </ul>
<b>Class of Antibody</b>	IgM, Monoclonals
<b>Source</b>	Human cell line
<b>Specificity</b>	<ul style="list-style-type: none"> <li>• ANTI-c 100% to c antigen</li> </ul>
<b>Standardisation</b>	Follows AABB and FDA guidelines
<b>Storage / Stability</b>	2-8°C, 24 months

Presentation		Pack	Cat.No.
Rh genotyping reagents	ANTI-c	2 ml	10171002
		5 ml	10171005