



 <b>ANTI-C3d</b>	<b>MONOSPECIFIC COOMBS SERA</b>
<b>Application</b>	Monospecific Coombs reagent for DAT and IAT
<b>Class of Antibody</b>	IgM, Monoclonal antibodies
<b>Specificity</b>	Specific to complement component C3d
<b>Standardisation</b>	Follows AABB and FDA guidelines
<b>Storage / Stability</b>	2-8°C, 24 months

<b>Presentation</b>	<b>Pack</b>	<b>Cat.No.</b>
<b>ANTI-C3d MONOSPECIFIC COOMBS SERA</b>	2 ml	10181002
	5 ml	10181005
	10 ml	10181010

