



ERYBANK[®]
ANTI-HUMAN IgG

MONOSPECIFIC COOMBS SERA

Application	Monospecific Coombs reagent for DAT and IAT
Class of Antibody	Polyclonal Anti-human IgG
Specificity	Specific to human gamma globulins
Standardisation	Follows AABB and FDA guidelines
Storage / Stability	2-8°C, 24 months

Presentation	Pack	Cat.No.
ANTI-HUMAN IgG MONOSPECIFIC COOMBS SERA	2 ml	10240002

