



ERYGEN-ABO	Red Blood Cell Panel for reverse grouping
Application	3% Reagent red blood cell panel for confirmation of the results of forward grouping
Specificity	ABO isoagglutinins
Reagent	Erygen-ABO panel of 3% reagent red blood cells comprises of: <ul style="list-style-type: none"> • Red blood cells of group A1 • Red blood cells of group B • Red blood cells of group O
Standardisation	Follows FDA guidelines
Storage / Stability	2-8°C, 50 days

Presentation	Pack	Cat.No.
Erygen – ABO	3x2ml	101950032

