



ERYGEN-AC

Coombs control cells

Application	3% Reagent red blood cell panel for quality control of anti-human globulin procedures
Source	Human
Specificity	Standardized to provide positive agglutination control for antihuman globulins containing IgG
Reagent	Erygen-AC panel of 3% reagent red blood cells comprises of: <ul style="list-style-type: none"> • Human erythrocytes of group O Rho (D) positive cells sensitized with Anti-D (IgG) monoclonal
Standardisation	Follows AABB and FDA guidelines
Storage / Stability	2-8°C, 50 days

Presentation	Pack	Cat.No.
Erygen – AC	2ml	10194302