



	<h2 style="text-align: center;">ABO Subgrouping Card Anti-A1</h2>
Application	For confirmation of A1 subgroup of A
Card Configuration	6 x Anti-A1 / Six profiles per card
Reagent	Matrix™ ABO Subgrouping Card Anti-A1 contains: <ul style="list-style-type: none"> • Anti-A1 lectin
Standardisation	Follows FDA guidelines
Storage / Stability	4-25°C, 12 months

Presentation	Pack	Cat.No.
Matrix™ ABO Subgrouping Card Anti-A1	24 Tests	102750024

