



	<h2 style="text-align: center;">ABO Subgrouping Card Anti-H</h2>
Application	For detection of H antigen
Card Configuration	6 x Anti-H / Six profiles per card
Reagent	Matrix™ ABO Subgrouping Card Anti-H contains: <ul style="list-style-type: none"> • Anti-H lectin
Standardisation	Follows FDA guidelines
Storage / Stability	4-25°C, 12 months

Presentation	Pack	Cat.No.
Matrix™ ABO Subgrouping Card Anti-H	24 Tests	102760024

