



	AHG (Coombs) Test Card
Application	For Direct Coombs test (DAT) and indirect Coombs test (IAT) including antibody screening, antibody identification and compatibility testing.
Card Configuration	6 x AHG / Profiles per card depends on the test parameter
Reagent	Matrix™ Anti-Human Globulin card contains: <ul style="list-style-type: none"> • Six microtubes prefilled with gel in suitable buffer containing polyspecific Coombs reagent (Anti- IgG + Monoclonal Anti-C3d)
Standardisation	Follows FDA guidelines
Storage / Stability	4-25°C, 12 months

Presentation	Pack	Cat.No.
Matrix™ AHG (Coombs) Test card	24 cards	10256024

