



	<h3>Coombs and Neutral Gel Card</h3>
Application	Combination of Coombs and neutral gel can be used for multiple applications such as antibody screening in AHG and Enzyme phase or Compatibility testing in AHG and Enzyme/saline phase
Card Configuration	3 x AHG-3 x Neutral / Profiles per card depends on test parameter
Reagent	Matrix™ Coombs and Neutral Gel Card contains: <ul style="list-style-type: none"> • Polyspecific Coombs (Anti-IgG+ monoclonal C3d) and Neutral gel
Standardisation	Follows FDA guidelines
Storage / Stability	4-25°C, 12 months

Presentation	Pack	Cat.No.
Matrix™ Coombs and Neutral Gel Card	24 Tests	102640024

