



	<h3>Coombs Anti-IgG Card</h3>
Application	For all Coombs test where anti complement activity is not required, like direct Coombs test and indirect Coombs test including antibody screening, antibody identification, compatibility testing and weak D testing
Card Configuration	6 x Anti-IgG / Profile per card depends on test parameter
Reagent	Matrix™ Coombs Anti-IgG Card contains: <ul style="list-style-type: none"> • Monospecific Coombs (Anti-IgG)
Standardisation	Follows FDA guidelines
Storage / Stability	4-25°C, 12 months

Presentation	Pack	Cat.No.
Matrix™ Coombs Anti-IgG Card	24 Tests	102660024

