



	<h3>Coombs Anti-IgG Card</h3>
<b>Application</b>	For all Coombs test where anti complement activity is not required, like direct Coombs test and indirect Coombs test including antibody screening, antibody identification, compatibility testing and weak D testing
<b>Card Configuration</b>	6 x Anti-IgG / Profile per card depends on test parameter
<b>Reagent</b>	Matrix™ Coombs Anti-IgG Card contains: <ul style="list-style-type: none"> <li>• Monospecific Coombs (Anti-IgG)</li> </ul>
<b>Standardisation</b>	Follows FDA guidelines
<b>Storage / Stability</b>	4-25°C, 12 months

Presentation	Pack	Cat.No.
<b>Matrix™ Coombs Anti-IgG Card</b>	24 Tests	102220024
	48 Tests	102220048

