



	<b>Diluent – 2 LISS</b>
<b>Application</b>	For preparation of red blood cell suspensions used for Matrix™ gel system
<b>Reagent</b>	Buffered low ionic strength solution of appropriate sodium chloride molarity
<b>Standardisation</b>	Follows FDA guidelines
<b>Storage / Stability</b>	2-8°C, 12 months

Presentation	Pack	Cat.No.
<b>Matrix™ Diluent-2 LISS</b>	100 ml	102570100
	250 ml	10257250
	500 ml	102570500

