



Direct Antiglobulin test card

Application	For detection and differentiation of IgG and or C3d activity
Card Configuration	IgG-C3d-Ctrl- IgG-C3d-Ctrl / two profile per card
Reagent	Matrix™ Direct Antiglobulin test Card contains: <ul style="list-style-type: none"> • Anti-Human IgG, Monoclonal Anti-C3d and neutral gel in appropriate microtubes.
Standardisation	Follows FDA guidelines
Storage / Stability	4-25°C, 12 months

Presentation	Pack	Cat.No.
Matrix™ Direct Antiglobulin test Card	24 Tests	102670024