




|   |   |
|---|---|
|  | <b>ERYGEN-ABO</b>   |
| <b>Application</b>  | 0.8% Reagent red blood cell panel for confirmation of the results of forward grouping   |
| <b>Specificity</b>  | ABO isoagglutinins  |
| <b>Reagent</b>  | Matrix Erygen-ABO panel of 0.8% reagent red blood cells comprises of: <ul style="list-style-type: none"> <li>• Red blood cells of group A1</li> <li>• Red blood cells of group B</li> <li>• Red blood cells of group O</li> </ul> |
| <b>Standardisation</b>  | Follows FDA guidelines  |
| <b>Storage / Stability</b>  | 2-8°C, 36 days  |

| Presentation       | Pack     | Cat.No.   |
|--------------------|----------|-----------|
| Matrix™ ERYGEN-ABO | 3 x 2 ml | 102070032 |

