



	ERYGEN-AS
Application	0.8% reagent red blood cells for Antibody Screening on Matrix™ Gel Cards
Reagent	Matrix™ ERYGEN-AS contains: <ul style="list-style-type: none"> • 3 cell panel of 0.8% reagent red blood cells for Antibody Screening on Matrix™ Gel Cards
Standardisation	Follows FDA guidelines
Storage / Stability	2-8°C, 35 days

Presentation	Pack	Cat.No.
Matrix™ ERYGEN-AS	3 x 2 ml	102700032

