



	ABO and Rho (D) Forward and Reverse grouping card with Autocontrol
Application	For determination of ABO and Rho (D) blood group systems and reverse grouping for confirmation of ABO grouping
Card Configuration	A-B-D-Ctrl-A1-B / One profile per card
Reagent	Matrix™ ABO and Rho (D) Forward and reverse grouping card with autocontrol contains: <ul style="list-style-type: none"> • Monoclonal Anti-A, Anti-B and Anti-D (VI-) prefilled with gel from microtube 1-3 • Microtube 4 contains neutral gel which serves as a negative control • Microtube 5 and 6 contains neutral gel for reverse grouping
Standardisation	Follows FDA guidelines
Storage / Stability	4-25°C, 12 months

Presentation	Pack	Cat.No.
Matrix™ ABO and Rho (D) Forward and reverse grouping card with autocontrol	24 Tests	10260024

