



	ABO and Rho (D) Forward Grouping & Crossmatch card
Application	For confirmation of ABO and Rho (D) blood group systems for both patient and donor, ABO and Rho (D) compatibility testing, major crossmatch in Coombs and enzyme phase along with Coombs phase autocontrol for the patient
Card Configuration	A-B-D-ENZ-AHG-AHG / One profile per card
Reagent	Matrix™ ABO & Rho (D) forward grouping & Crossmatch card contains: <ul style="list-style-type: none"> • Monoclonal Anti-A, Anti-B and Anti-D (VI-) prefilled with gel from microtube 1-3 • Microtube 4 contains neutral gel for Enzyme phase major crossmatch • Microtube 5 and 6 contains polyspecific Coombs (Anti-IgG+ monoclonal C3d)
Standardisation	Follows FDA guidelines
Storage / Stability	4-25°C, 12 months

Presentation	Pack	Cat.No.
Matrix™ ABO & Rho (D) Forward Grouping & Crossmatch Card	24 Tests	102650024

