



	<h2>Forward Grouping DAT Card</h2>
Application	For detection of ABO and Rho (D) blood group systems using two different clones of Anti-D in separate microtubes and Direct Coombs test (DAT)
Card Configuration	A-B-D-D-ctrl-AHG / One profile per card
Reagent	Matrix™ Forward Grouping DAT Card contains: <ul style="list-style-type: none"> • Monoclonal Anti-A, Anti-B, Anti-D (VI-) and polyspecific Coombs (Anti-IgG + monoclonal C3d)
Standardisation	Follows FDA guidelines
Storage / Stability	4-25°C, 12 months

Presentation	Pack	Cat.No.
Matrix™ Forward Grouping DAT Card	24 Tests	102730024

