



	ABO/Rho (D)/AHG Neonate Group card
Application	For determination of ABO and Rho (D) status and direct Coombs test of neonates
Card Configuration	A-B-AB-D-Ctrl-AHG / One profile per card
Reagent	Matrix ABO/ Rho (D)/AHG Neonate Group card contains: <ul style="list-style-type: none"> • Six microtubes prefilled with gel in suitable buffer containing specific monoclonal antibodies Anti-A, Anti-B, Anti AB, Anti-D (VI-) and polyspecific Coombs (Anti-IgG + monoclonal C3d) and neutral gel (negative control) in appropriate positions
Standardisation	Follows FDA guidelines
Storage / Stability	4-25°C, 12 months

Presentation	Pack	Cat.No.
Matrix™ ABO/Rho (D)/AHG Neonate Group card	24 Tests	10259024

