



	<h2>Neutral Gel Card</h2>
Application	For saline and Enzyme phase tests including antibody screening, antibody identification, compatibility testing and ABO reverse grouping
Card Configuration	6 x Neutral gel / Profile per card depends on test parameter
Reagent	Matrix™ ABO reverse grouping card contains: <ul style="list-style-type: none"> • Neutral gel in all six microtubes
Standardisation	Follows FDA guidelines
Storage / Stability	4-25°C, 12 months

Presentation	Pack	Cat.No.
Matrix™ Neutral Gel card	24 cards	102740024

