



 Forward and Reverse Grouping Card with Sub Grouping	Matrix™ Octoplus Forward and Reverse Grouping Card with Sub Grouping
Application	For determination of ABO and Rho (D) blood group systems with subgrouping and reverse grouping for confirmation of ABO grouping
Card Configuration	A-B-D-Ctrl-A1-H-A1-B / One profile per card
Reagent	Matrix™ Octoplus Forward and Reverse Grouping Card with Sub Grouping contains: <ul style="list-style-type: none"> • Monoclonal Anti-A, Anti-B, Anti-AB and Anti-D (VI-) prefilled with gel from microtube 1-3 • Microtube 4 contains neutral gel which serves as a negative control • Microtube 5 contains Anti-A1 lectin and Microtube 6 contains Anti-H lectin • Microtube 7 and 8 contains neutral gel for reverse grouping
Standardisation	Follows FDA guidelines
Storage / Stability	4-25°C, 12 months

Presentation	Pack	Cat.No.
Matrix™ Octoplus Forward and Reverse Grouping Card with Sub Grouping	24 Cards	102560192

