




|  |  |
|--|--|
|  | <b>QAS level - II</b>  |
| <b>Application</b>   | Standardized, lyophilized monoclonal Anti-D IgG in specified concentration for Quality Assurance of reagents and techniques related to Matrix™ Gel System and classical tube technique |
| <b>Reagent</b>   | <b>Matrix™ QAS level -II</b> contains:<br>Standardized Anti-D IgG (Lyophilized)  |
| <b>Standardisation</b>   | Follows FDA guidelines   |
| <b>Storage / Stability</b>   | Lyophilized: 2-8°C, 12 months<br>Reconstituted: 2-8°C, one week  |

| Presentation          | Pack                              | Cat.No.   |
|-----------------------|-----------------------------------|-----------|
| Matrix™ QAS level -II | 0.5 ml Anti-D and<br>10 ml 6% BSA | 102780002 |

