


| | |
|--|--|
|  | <h2 style="text-align: center;">Rh Phenotype Card with Anti-D</h2> |
| Application | For determination of antigens of Rh system |
| Card Configuration | C-c-D-E-e-Ctrl / One Profile per card |
| Reagent | Matrix™ Rh Phenotype Card with Anti-D contains: <ul style="list-style-type: none"> • Monoclonal Anti-C, Anti-c, Anti-D (VI-), Anti-E, Anti-e with Neutral gel |
| Standardisation | Follows FDA guidelines |
| Storage / Stability | 4-25°C, 12 months |

| Presentation | Pack | Cat.No. |
|--|----------|-----------|
| Matrix™ Rh Phenotype Card with Anti-D | 24 Tests | 102630024 |

