



	<h2 style="text-align: center;">Rh Phenotype Card with Anti-D</h2>
Application	For determination of antigens of Rh system
Card Configuration	C-c-D-E-e-Ctrl / One Profile per card
Reagent	Matrix™ Rh Phenotype Card with Anti-D contains: <ul style="list-style-type: none"> • Monoclonal Anti-C, Anti-c, Anti-D (VI-), Anti-E, Anti-e with Neutral gel
Standardisation	Follows FDA guidelines
Storage / Stability	4-25°C, 12 months

Presentation	Pack	Cat.No.
Matrix™ Rh Phenotype Card with Anti-D	24 Tests	102200024
	48 Tests	102200048

