



Rh Phenotype Card with Anti-K

Application	For determination of Rh Phenotypes along with K antigen of Kell system
Card Configuration	C-c-E-e-K-Ctrl / One Profile per card
Reagent	Matrix™ Rh Phenotype Card with Anti-K contains: Monoclonal Anti-C, Anti-c, Anti-E, Anti-e, Anti-K and neutral gel in appropriate microtubes.
Standardisation	Follows FDA guidelines
Storage / Stability	4-25°C, 12 months

Presentation	Pack	Cat.No.
Matrix™ Rh Phenotype Card with Anti-K	12 Tests	102190012
	24 Tests	102190024

