



RHOFINAL[®]

ANTI-D

Application	Reagent for routine typing and confirmation of Rh ₀ (D) phenotypes by slide and tube test
Titre	≥ 1:256
Class of Antibody	IgM + IgG, Monoclonal Blend
Source	Human B cell line, EBV transformed
Specificity	100% to Rh ₀ (D) antigen
Standardisation	Follows AABB and FDA guidelines
Storage / Stability	2-8°C, 24 months

Presentation	Pack	Cat.No.
Rh₀(D) typing reagent	5 ml	101660005
	10 x 5 ml	101660105
	10 ml	101660010
	6 x 10 ml	101660010

