



	<h2>ReaCell Panel</h2>
Application	0.8% reagent red blood cells for Antibody Screening on Matrix™ Gel Cards
Reagent	ReaCell Panel contains: <ul style="list-style-type: none"> • 11 cell panel of 0.8% reagent red blood cells for Antibody Screening on Matrix™ Gel Cards
Standardisation	Follows FDA guidelines
Storage / Stability	2-8°C, 53 days

Presentation	Pack	Cat.No.
ReaCell Panel	11 x 5 ml	IOTH-0241

