



Calcium Chloride

Recalcification reagent for use with APTT reagents

Intended Use: _____

Calcium Chloride, a ready to use 0.025 M recalcification solution for APTT test.

Calcium Chloride kit components:

Reagent	Ready to use Liquid, 0.025 mol/L CaCl ₂
Other Accessories	Package Insert

FEATURES	BENEFITS
Ready to use	Ready to use Liquid, 0.025 mol/L for APTT test
Convenient pack sizes	Ideal for all users
Precisely prepared reagent	Ensures that optimum results are obtained in each case
Easy to use	No measuring hassles

Storage / Stability	Temperature	Duration
Unopened kit	20-25°C	24 Months
Opened kit	20-25°C	24 Months

Available Pack Sizes	
10 ml	10 x 10 ml