



Factor VIII Deficient Plasma

Intended Use: _____

For use as a deficient substrate in coagulation assays

Factor VIII Deficient plasma Kit components:

Reagent	Freeze dried, stabilized Factor VIII immune depleted plasma with Factor VIII activity < 1 % & other coagulation factor activity in the normal range
Other Accessories	Package Insert

FEATURES	BENEFITS
Freeze Dried Reagent	Aids in minimizing microbial contamination
Stable reagent	Freeze dried, stabilized citrated plasma
Precisely Prepared Reagent	When reconstituted appropriately, ensures precise results
Convenient Pack size	Ideal for all users

Storage / Stability	Temperature	Duration
Unopened kit	2-8°C	18 months
Reconstituted reagent	2-8°C	8 hours
	-20°C	1 month

Available Pack Size
3 x 1 ml