



Plasmatrol H - I/II/III
Control Plasma for coagulation assays

Intended Use: _____

Plasmatrol H - I/II/III, control plasma for calibration and quality control of PT, APTT, TT, Fibrinogen.

Plasmatrol H - I/II/III kit components:

Reagent	Freeze dried, stabilized citrated plasma
Other Accessories	Package Insert, Assay value sheet

FEATURES	BENEFITS
Stable reagent	Freeze dried, stabilized citrated plasma
Convenient pack size	Ideal for all users
Precisely prepared reagent	When diluted appropriately, ensures precise results
Freeze dried reagent	Aids in minimizing microbial contamination

Storage / Stability	Temperature	Duration
Unopened kit	2-8°C	18 months
Reconstituted reagent	2-8°C	8 hours
	25-30°C	3 hours

Available Pack Sizes
6 x 1 ml