



Plasmatrol R
Reference Plasma for use with coagulation assays

Intended Use: _____

Plasmatrol R, Reference Plasma for calibration and quality control of PT, APTT, Fibrinogen, Factors II, V, VII, VIII, IX, X, XI, XII.

Plasmatrol R kit components:

Table with 2 columns: Component, Description. Rows: Reagent (Freeze dried, stabilized citrated plasma), Other Accessories (Package Insert, Assay value sheet)

Table with 2 columns: FEATURES, BENEFITS. Rows: Stable reagent, Convenient pack size, Precisely prepared reagent, Freeze dried reagent

Table with 3 columns: Storage / Stability, Temperature, Duration. Rows: Unopened kit, Reconstituted stability (25-30°C, 2-8°C)

Table with 1 column: Available Pack Sizes. Row: 6 x 1 ml