



Profact 3.2% Buffered Tri-Sodium Citrate

Intended Use: _____

Profact, a ready to use 3.2% Buffered Tri-Sodium Citrate solution for coagulation assays and ESR by Westergren method.

Profact kit components:

Reagent	3.2% Buffered Tri-Sodium Citrate
Other Accessories	Package Insert

FEATURES	BENEFITS
Convenient pack sizes	Ideal for all users
Precisely prepared reagent	Ensures that optimum results are obtained in each case
Easy to use	No measuring hassles
Incorporates 3.2% sodium citrate in an unique protective solution	Arrests shift in pH of plasma during centrifugation

Storage/Stability	Temperature	Duration
Unopened kit	2-8°C	18 Months
Opened kit	2-8°C	18 Months

Available Pack Sizes
20 ml