



	<b>For use as a deficient substrate in coagulation assays</b>
<b>Application</b>	<p>A Particular Factor deficient plasma is useful in the quantitative determination of activity of that factor using appropriate clot based assay. The activity of the factor VII of Extrinsic pathway is performed by the PT test.</p>
<b>Reagent</b>	<p>Factor VII immune depleted plasma with Factor VII activity &lt; 1 % &amp; other coagulation factor activity in the normal range</p>
<b>Storage / Stability</b>	<p>2-8°C for 18 months</p>
<b>In Use Stability</b>	<p>Once reconstituted factor deficient plasma is stable for</p> <ul style="list-style-type: none"> <li>• 8 hours when stored at 2-8°C</li> <li>• 1 month at -20°C if rapidly frozen</li> </ul>

Presentation	Pack	Cat.No.
Factor VII deficient Plasma	3 x1 ml	110560031