



## Thrombin time Reagent

<b>Application</b>	For the qualitative estimation and evaluation of Fibrinogen in plasma
<b>Reagent</b>	Lyophilised Thrombin reagent with appropriate preservative
<b>Source</b>	Bovine
<b>Normal Range</b>	5-15 seconds
<b>Storage / Stability</b>	2-8°C, 18 months
<b>In Use Stability</b>	Reconstituted, <ul style="list-style-type: none"> <li>• 2-8°C, 6 days</li> <li>• 20-25°C, 4 hours</li> </ul>

Presentation	Pack	Cat.No.
FIBROSCREEN	2 ml	10640002
	6 x 1 ml	10640061