



APTT Reagent with 3.2% Buffered Citrate blood collection anticoagulant

Application	Combipack of Liquicelin-E (APTT) reagent and 3.2% citrate with zwitterionic buffer for achieving precise, consistent and accurate results in APTT determination.
Reagent	<ul style="list-style-type: none"> • Ready to use Liquid stable, Cephaloplastin reagent with Ellagic acid as an activator • 3.2% Buffered Tri-sodium citrate
Source	Rabbit Brain Phospholipids
Normal Range	22-34 seconds
Storage / Stability	2-8°C, 18 months

Presentation	Pack	Cat.No.
LIQUICELIN-E® PRECISION COMBIPACK contains Liquicelin-E® and 3.2% buffered citrate	2.5 ml (25 Tests)	10631025

