



	<p>Thromboplastin Reagent with 3.2% Buffered Citrate Blood collection anticoagulant</p>
<p>Application</p>	<p>Combipack of Liquiplastin (PT) reagent and 3.2% citrate with zwitterionic buffer for achieving precise, consistent and accurate results in PT determinations</p>
<p>Reagent</p>	<ul style="list-style-type: none"> • Ready to use Liquid stable, Calcified Thromboplastin reagent • 3.2% Buffered Tri-sodium citrate
<p>Source</p>	<p>Rabbit Brain</p>
<p>Normal Range</p>	<p>10 - 14 seconds</p>
<p>ISI value</p>	<p>Optimised at ~ 1.5 (1.5 -1.7)</p>
<p>Storage / Stability</p>	<p>2-8°C, 18 months 18-25°C, 1 week 37°C, 2 days</p>

Presentation	Pack	Cat.No.
<p>LIQUIPLASTIN® PRECISION COMBIPACK contains Liquiplastin and 3.2% buffered citrate</p>	<p>5 ml (25 Tests)</p>	<p>10611025</p>

