



## Thromboplastin Reagent

<b>Application</b>	For the Prothrombin Time determination in plasma, monitoring Oral Anticoagulant therapy and assessment of Extrinsic Pathway
<b>Reagent</b>	Lyophilized Calcium Thromboplastin reagent
<b>Source</b>	Rabbit Brain
<b>Normal Range</b>	11-15 seconds
<b>ISI value</b>	Optimised at ~ 1.0
<b>Storage / Stability</b>	2-8°C, 18 months
<b>In Use Stability</b>	Reconstituted stability, 10 days at 2-8°C

Presentation	Pack	Cat.No.
LYOPLASTIN®	12 x 4 ml	10623124