



	<h2>Control Plasma</h2>
Application	Normal Control Plasma for calibration and Quality Control of PT, APTT, TT, Fibrinogen.
Reagent	Freeze dried, stabilized citrated plasma
Storage / Stability	2-8°C, 12 months
In Use Stability	Once reconstituted should be used within 3 hours <ul style="list-style-type: none"> • 8 hours when stored at 2-8°C • 3 hours when stored at 25-30°C

Presentation	Pack	Cat.No.
PLASMATROL-H I	6 x 1 ml	11040061

