



	<b>Control Plasma</b>
<b>Application</b>	Abnormal Control Plasma for calibration and Quality Control of PT, APTT, TT, Fibrinogen.
<b>Reagent</b>	Freeze dried, stabilized citrated plasma
<b>Storage / Stability</b>	2-8°C, 12 months
<b>In Use Stability</b>	Once reconstituted should be used within 3 hours <ul style="list-style-type: none"> <li>• 8 hours when stored at 2-8°C</li> <li>• 3 hours when stored at 25-30°C</li> </ul>

Presentation	Pack	Cat.No.
PLASMATROL-H III	6 x 1 ml	11042061

