



	<h2>Reference Plasma</h2>
<b>Application</b>	Reference Plasma for calibration and Quality Control of PT, APTT, Fibrinogen, Factors II, V, VII, VIII, IX, X, XI, XII.
<b>Reagent</b>	Freeze dried, stabilized citrated plasma
<b>Storage / Stability</b>	2-8°C, 12 months
<b>In Use Stability</b>	Once reconstituted should be used within 3 hours <ul style="list-style-type: none"> <li>• 7 hours when stored at 2-8°C</li> <li>• 3 hours when stored at 25-30°C</li> </ul>

Presentation	Pack	Cat.No.
PLASMATROL-R	6 x 1 ml	11043061

